

Ardagh Dental
Supply

Technical
Solutions

**SPRING/SUMMER
PROMOTIONS
2016**

bredent

Thinking of flexible dentures??

visio.lign

Clear shapes &
absolute stability



Training Available



More
quality of life
with metal free
restorations

Details
Available



NEW

IPS Ivocolor

Universal stains and glazes

One for all



TASSO SAMBANIDIS 0872847795
info@ardaghdental.ie

ARDAGH DENTAL (01)8551600

ALL PRICES EXCLUDE VAT & OE

Crown & Bridge

IPS Style®

The first metal-ceramic containing patented oxyapatite

**MAKE IT
STYLE YOUR!**

Make it fast. IPS Style's low shrinkage saves you time.

Make it easy. The material's ease of use ensures reliable results.

Make it natural. The patented oxyapatite controls the translucency and depth effect.

IPS Ivocolor® – One for All

Universal shades and stains

IPS Ivocolor Shade

Shade adaptation: swift, easy, exquisite

Shade Pastes - 3g Refill

€29.00



IPS Ivocolor Shade Dentine



Shade Incisal



667651 Shade Kit

€495.00

IPS Ivocolor Essence

Effects, nuances and accentuations:
faithfully reproducing nature

Essence Powders - 1.8g Refill

€12.50

Glaze Pastes - 3g Refill

€29.00

Glaze Powders - 1.8g Refill

€9.95



667650 Essence Kit

€395.00

IPS Ivocolor Mixing Liquid Allround

€24.00

IPS Ivocolor Mixing Liquid Longlife

€24.00

IPS Ivocolor Essence

Essence Basic



NEW PRODUCT

IPS Ivocolor: Glaze

The new standard in glazing materials

IPS Ivocolor Glaze for reliable firing results with exact shades:

- High gloss in a single firing process
- Consistent gloss, independent of the ceramic substrate (e.g. cut-back)
- No grey or white discolouration
- Glazes with and without fluorescent effect
- Available in a powder and paste delivery form



**Ask your Product Specialist
for all exchange offers**


ORDER NOW

IPS Style®

The first metal ceramic containing patented oxyapatite

NEW

MAKE IT YOUR STYLE!



"Too easy!"

Ask your Product Specialist
for all exchange offers

ORDER NOW



Ceram

Press

One

IPS Style® Ceram Starter Kit

IPS Style® Ceram A2 Intro Kit

*Everything you need for
an A2 restoration!*

€795.00



€140.00



**ivoclar
vivadent®**



IPS e.max Press

Multi ingots - 5 pcs €95.00

HT ingots - 5 pcs €77.00

MT ingots - 5 pcs €77.00

MO ingots - 5 pcs €77.00

HO ingots - 5 pcs €77.00

Impulse ingots - 5 pcs €77.00



IPS e.max® Press Multi - Intro Kit A-D

Kit contains:

- 5 pcs IPS e.max Press Multi A1
- 5 pcs IPS e.max Press Multi A2
- 1x IPS Multi Investment Ring Base 200g
- 80 pcs IPS Multi Wax Pattern Form A
- 80 pcs IPS Multi Wax Pattern Form B
- 1x IPS Multi Sprue Guide 200g
- 10 pcs IPS Multi One-Way-Plunger



★ **€318.00**

IPS e.max® Press Multi - One-way Plunger



★ **€150.00**

IPS e.max® Press Multi - Ingot Refills

- 10 Shades (1 Bleach BL, 9 A-D shades)
BL2, A1, A2, A3, A3.5, B1, B2, C1, C2, D2



★ **€95.00**

SPECIAL OFFER

**IPS e.max® Press Multi
1 day Workshop**
ASK FOR DETAILS

IPS e.max® Press Multi - Multi Wax Patterns

- Pattern Form A (green) or
- Pattern Form B (pink)



★ **€35.50**

IPS e.max® Press Multi - Investment Ring System 200g

- Base & Gauge



★ **€33.50**

IPS e.max® Press Multi - Sprue Guide 200g

- IPS Multi Investment Ring Base not included

★ **€13.00**



Furnaces

Ivoclar Vivadent has a range of innovative ceramic and press furnaces that are sure to assist in the productivity of your laboratory. Features like*:

- Power Saving Technology
- Automatic double-range temperature calibration
- Integrated infrared thermal camera
- Multimedia functions

* Not on all models

Programat® EP 3010

The user-friendly press and ceramic furnace



€6,695 Special Offer Request

Programat® EP 5010

The intelligent press and ceramic furnace

€8,995 Special Offer Request

Programat® P310

The optimal furnace for the daily laboratory routine



€3,395 Special Offer Request

Programat® P510

Excellent results, with infrared camera



€4,495 Special Offer Request

Programat P710®

Packed with innovative features, the new Programat P710 can do more than just fire.



€5,495 Special Offer Request

* Discount applied when booking and paying. No cash alternative.

Prices exclude GST, Freight & Installation, and may change without notice. E & OE

Boiling Out Units

Dreve

€3800

£3300

FREE

Package A contains

- Labormat SD
- 4 flask baskets
- 1 tub Unisol E (2.5 kg)
- 3 canisters Unisol W (5 l)
- 10 filter mats

FREE

Package B contains

- Labormat TH
- 2 flask baskets
- 2 tubs Unisol E (2.5 kg)
- 4 canisters Unisol W (5 l)
- 10 filter mats

Dreve

Hot and cold solid Pressure polymerisation unit

€1695.00

Made for Composites.

Dreve

€1250.00

Dreve

€1595,00

Dosper evo

Automatic mixing and dosing unit for thin-flowing, addition-vulcanising duplication silicones. The working process becomes shorter and safer. Materials reach optimum quality values. The unit is constantly ready for use and very economical, a profitable machine also for smaller lab.

Lightplast base plates

Dreve

Light-curing base plates

Pink / Natural

€51.95 (50 pcs)

Dreve

MultiSpot®

Mobile curing unit for dental technicians

€259,00

**SPECIAL
OFFER**

Dreve

StoneBite - Bite Registration Material

€35.95

2 x 50ml
Double Cartridges
incl. 12 mixing tips

Dreve

NanoVarnish

Light-curing one-component
lacquer on the basis of
especially structured nano-
particles.

€49.95
(20ml)



NanoVarnish
Licht-aktiver NanoVarnish
nano-lackier photosensibilisiert
Photopolymersubstanz
Contains: Methacrylate
Polyethylene
Content: 20 ml
Item no.: 18905

Drewe Dentamid GmbH
Max-Ramke-Str. 91 · 50623 Uckerath
Tel./Phone: +49 2303 8807-0

**ivoclar
vivadent**
passion vision innovation



PRODUCT HIGHLIGHT

Teeth Selection

Tooth Consignments Available Now
Please speak to your local Product Specialist



Take advantage of this fantastic offer to receive 10, 20 or 50 sets of teeth free of charge*.

This offer is available on all Ivoclar Vivadent tooth ranges, teeth purchased can only be from one brand: (50 x SR Phonares® II)



**ivoclar
vivadent**
passion vision innovation



PHONARES 2

POSTERIOR € 25.00
ANTERIOR € 32.00



**ivoclar
vivadent**
passion vision innovation



VIVADENT PE

POSTERIOR € 13.20
ANTERIOR € 18.20



**ivoclar
vivadent**
passion vision innovation




IVOSTAR / GNATHOSTAR

POSTERIOR € 5.90
ANTERIOR € 5.90



**ivoclar
vivadent**
passion vision innovation



VIVODENT DCL

POSTERIOR € 16.00
ANTERIOR € 20.00

TASSO SAMBANIDIS 0872847795
info@ardagh dental.ie

ARDAGH DENTAL (01)8551600
ALL PRICES EXCLUDE VAT & OE


SR Vivodent® S PE

Inspired by nature, prosthetically optimised




Features	Customer Benefits
Design & texture	<ul style="list-style-type: none"> New and modern moulds Circular enamel appearance <p>For an extraordinary mould design</p>
Brilliant shades	<ul style="list-style-type: none"> Exclusive PE layering High shade intensity <p>For a unique life-like effect</p>
Material	<ul style="list-style-type: none"> Highly cross-linked DCL material Higher compressive strength compared to PMMA and similar flexibility <p>Wide range of indications</p>


SMALL



MEDIUM



LARGE



Tooth Consignment
with your first tooth consignment of SR Vivodent S PE teeth you can receive up to 150 sets
FREE

Set of Anteriors €18.50
Set of Posteriors €13.40
Sample Set €26.25


The Multi-functional Shade Guide

Integrated facialmeter

- Integrated into the holder of the shade guide
- Used to determine the interalar width.
The measured value is then assigned to a size (small, medium or large) on the flip side of the facial meter

Integrated papillameter

- Integrated in the flip side of the shade guide in a contrasting colour
- Used to determine the length of the lip closure line. It is the first reference for the occlusal plane and the length of the anterior teeth
- Autoclavable



Request a copy of the new Technical Documentation TODAY

IvoBase® System

The innovative denture base



SPECIAL OFFER AVAILABLE

IvoBase® Hybrid material is characterised by its short processing times. Within 35 minutes, the material thoroughly cures in the IvoBase® Injector. IvoBase® High Impact material is characterised for fracture toughness and can be polymerised in 50 minutes.



OptraSculpt® Pad

OptraSculpt Pad is a contouring instrument with special foam pad attachments, designed for the efficient, non-stick forming and shaping of composites.

- Professional aesthetic results
- Creation of smooth and even surfaces

ASSORTMENT - €79.00

REFILL - €33.00

**BUY 1 GET 1 FREE
FIRST TIME USER OFFER**



SR Nexco Paste

Lab composite for a lifelike appearance

SR Nexco Paste, the purely light-curing lab composite, embodies aesthetics, shade tolerance and versatility. The composite contains micro-opal fillers, which provide natural-looking opalescence and translucency.

SR Nexco restorations exhibit a lifelike appearance due to the micro-opal fillers. In addition, the restorations feature high shade stability and are therefore durable.

In combination with the shading that is tolerant of varying layer thicknesses, true-to-nature optical properties can be achieved for fixed and removable dentures also involving artificial gingiva.



SR Nexco Paste Starter A-D

€480

SR Nexco Paste Gingiva Kit

€195

Refills €17.00



Removable Prosthetics

IvoBase® Hybrid

For all removable dentures

- High accuracy of fit due to shrinkage compensation
- Outstanding bond with resin teeth
- Polymerisation time: only 35 minutes
- **Shades:** Pink; Pink-V; Preference; Clear
With heightened opacity: Pink-V Implant; Preference Implant

20x Capsules

€149.00



IvoBase® High Impact

For immediate dentures & implant retained dentures

- Highly impact-resistant and fracture-proof material
- High accuracy of fit due to shrinkage compensation
- Polymerisation time: 50 minutes
- For exceptionally high stability requirements
- **Shades:** Pink; Pink-V; Preference
With heightened opacity: Pink-V Implant; Preference Implant; 34-V

20x Capsules

€169.00



SR Ivocap® High Impact

For complete dentures, partial dentures and relines

- Precision fit without vertical increase of occlusion
 - Outstanding bond with resin teeth
 - Excellent polishability
 - Tried-and-tested injection technique
- Shades: Pink, Pink-V, US-L, US-P, S, Preference, US-D, Pink-V Implant, Preference Implant

50x Capsules

€210.00



SR Ivocap® Clear

For orthodontic appliances, palatal plates and bite guard splints

- Excellent transparency
- Shade: Clear

50x Capsules

€195.00



SR Ivocap® Elastomer

For positioners, mouth guards, soft bite guards and snore guard applications

- Comfortable to wear
 - High elasticity and tear resistance
 - Neutral taste and smell
- Shade: Clear

10x Capsules

€125.00



ProBase® Cold

Suitable for both pouring and packing techniques. Its excellent properties ensure long-lasting denture bases. ProBase Cold is characterised by its outstanding flow and moulding properties. Ideal for repairs.

- Shades: US-P, Clear, Preference, P-V, P-V Implant

Standard Pack
2x 500g Powder/
500ml Monomer

€150.00



€86.00
(2 x 500g)

€50.00
(500ml)

ProBase® Hot

Suitable for partial and combination dentures (conventional or implant-supported), complete dentures, hybrid dentures as well as for relines of conventional and implant-supported dentures.

- Shades: US-P, US-L, Pink, P-V, Clear

Standard Pack
2x 500g Powder/
500ml Monomer

€85.00



Prosthetics



Centric Tray

- Reliable determination of the centric relation
- Subsequent adjustments are reduced
- Reusable, hence time-saving and economical
- For indications in both complete and partial denture prosthetics conventional or implant supported
- Autoclavable

€43.00



Stratos 100

- Average Value

- Fixed protrusive path of 30°
- Bennett angle of 15°
- Economically priced
- Fully integrated into the Stratos system
- 576820 (as pictured)

€317.00



Gnathometer

Gnathometer "M" (incl. UTS Adapter) is an intraoral registration device that enables you to determine the occlusal position of edentulous patients.

Needle-point tracing is used to check the occlusal position quickly and easily. The precision of this technique is particularly important in the fabrication of implant-supported dentures.

€315.00



Stratos 200

- Semi Adjustable

- Universally applicable articulator for both fixed and removable prosthesis
- Large model mounting space
- Easy open access to models
- Expandable system with a wide range of accessories
- 536381 (as pictured)

€595.00



Universal Transferbow System - 3D

€595.00

The UTS 3D transfer system is suitable for TMJ-related model transfer into articulators.

With the 3D registration joint, the registration data can be fixed with only one screw in just a few seconds. Due to the adjustable nosepiece, the UTS 3D can be universally set according to the Camper's Plane (CP) or the Frankfurt Horizontal (FH).



Stratos 300

- Fully Adjustable

- Designed to provide a full range of condylar adjustments
- Ability to personalise the settings of the articulator to better simulate the patient's jaw movements
- 576821 (as pictured)

€745.00



Personalised Assortment

- Plane indicator
- Bite fork support
- CE registration joint holder
- Incisal plates
- 590221

€255.00



Average Value Assortment

- Instrument carrier
- Horizontal guide
- 2D setting-up template
- Set-up table
- 536382

€450.00



3D setting-up template

For correct tooth set-up of personalised models in conjunction with the transferbow taking the Spee and Wilson curves into account.

- 536398

€319.00



TASSO SAMBANIDIS 0872847795

ARDAGH DENTAL (01)8551600

ALL PRICES EXCLUDE VAT & OE

Retention Sil

Complete Set

Designation	Shore hardness	Pull-off forces	REF
retention.sil 200	25 SH	200 g / 2 Newton	580 RTS 25
retention.sil 400	50 SH	400 g / 4 Newton	580 RTS 50
retention.sil 600	65 SH	600 g / 6 Newton	580 RTS 65



Multi-Sil-Primer
5 ml
REF 520 D100 4



Special Silicon Trimmer
Ø 4,1 mm
REF SKY-DR 41

The special silicone cutter is suitable for angled and lab handpiece.



REF 580 RT SET
Retention.sil in 3 hardnesses
Mixing Tips and application
instruments

€110



Dental laboratory sandblaster /
pneumatic / 2 tanks
PUNA

€1,295



Dental Manufacturing S.p.A.

We take care of your smile

ANTERIOR AND POSTERIOR

€3.50

BUY 6 SETS GET 1 SET FREE



*FREE SHADE GUIDE
AVAILABLE*

Reitel Steamy mini



€995

Acry RockV

Glass Beads

Blasting material for gentle cleaning and
compacting of sensitive surfaces of precious
and non-precious alloys
for example crowns, bridges and model casting.

Major features at a glance:

made in Germany
particular suitability: gentle cleaning and
compacting of sensitive surfaces
silicosis-free
available as corn class 50 µm and 125 µm



Aluminium oxyde

Blasting material white: very efficient
blasting caused by enormous hardness
and sharp-edged grains.

Major features at a glance:

made in Germany
purity grade 99.8 %
complies with industrial safety regulations
silicosis-free
3 sorting grades: 50 µm - 110 µm - 250 µm




50µm - 8kg €32
50µm - 15kg €52
125µm - 8kg €32
125µm - 15kg €52



50µm - 10kg €42
50µm - 20kg €73
110µm - 10kg €42
110µm - 20kg €73
250µm - 10kg €42
250µm - 20kg €73

bredent

thermopress 400 machine
€4,850



More
quality of life
with metal free
restorations



- Five biocompatible systems in one
- No compressed air
- Metal free
- Fully automatic
- Presses two flasks at once
- Denture frameworks clasps crowns & bridges



Accessories pack €765.00

Operation Equipment : Tools, Flask, Thermopaste, etc

Expando Rock €80.00

Complete System €5,730

Consumables

Ardagh | Dental
Supply
Co. Ltd.



Single cartridge

16g - 20g €14.50
24g €16.00
30g €17.50

Boxes of cartridges

20 x 16g €265.00
20 x 20g €285.00
20 x 24g €310.00
20 x 30g €335.00

Bulk buy (500g)

1 x 500g €140.00

TRAINING COURSE AVAILABLE

At the Irish Centre of Dental Education



visio.lign

Clear shapes & absolute stability



**COURSES
AVAILABLE**

buy with a free place on training course

- Veneering system for processing crowns, bridges and dentures
- 7 year stain resistance guarantee
- Ideal for working on implants, bars and chrome frameworks
- Full system includes veneers, teeth, customisation composites, opaques, gums and much more
- An alternate to composite build-ups

bre.Lux PowerUnit 2

€1,800



bre.Lux hand light

€475



Novo.lign two shade kit

€1,850

anterior & posterior, upper & lower, 2 shade, cabinet & all accessories

€2,400

Novo.lign/Neo.lign two shade kit

veneers & teeth, anterior & posterior, upper & lower, 2 shade, cabinet & all accessories

Novo.lign/Neo.lign two shade kit+


bre.lux2

veneers & teeth, anterior & posterior, upper & lower, 2 shade, cabinet & all accessories

€3,995

visio.lign ToolKit

€140.00



visio.lign

Erkoform-3d^{motion}

Erkoform-3d^{motion}

Vacuum thermoforming unit with touchless temperature control (EP 1 905 380), pre-vacuum, touchpanel and automated thermoforming process

€1795



Erkoform-3d +

Erkoform-3d +
Vacuum thermoforming unit with touchless temperature control (EP 1 905 380), pre-vacuum and touchpanel

€1550



SPECIAL OFFER
€150 worth of materials
with the purchase of any
Erkoform machine

OCCLUFORM-3



€450

Occludator integrated in a thermoforming unit for imprinting the opposing bite. The Occluform-3 can only be installed in Erkoform units (-3/-3d motion/-3d+)


**UPCOMING COURSE DATES FOR
HANDS-ON INTRODUCTION TO PLAY
SAFE SPORTS GUARDS, BLEACHING
TRAYS AND NIGHT SPLINTS**

**BOOK NOW
AVAILABLE DATES
11TH/12TH/13TH/14TH
OCTOBER 2016**

ERKOSKIN

Spacer for bleaching trays

€17.95
(50ml)




Finishing Set Quick 2

€59.95

Lisko-S



€25.95



Take Off Pliers



€35.95



Erkodent Hot Air Burner

€75.95




Ardagh Dental Supply

TASSO SAMBANIDIS 0872847795

ERKODENT®

ERKOFLEX BLEACH



€23.50 - 1mm (20pc)
 €93.95 - 1mm (100pc)
 €15.85 - 1.5mm (10pc)
 €21.85 - 2mm (10pc)

ERKOLOC / PRO



BITE SPLINTS

ERKOLOC

€32.25- 3mm (10pc)
 €128 - 3mm (50pc)

ERKOLOC pro soft/hard

€41.85- 3mm (10pc)
 €169.60- 3mm (50pc)
 €62.27- 1.3mm (20pc)

ERKODUR



**ALL FOILS AVAILABLE
IN SQUARE**

Play Safe Mouth Guards



RETAINERS

€19.62- 0.6mm (20pc)
 €29.80 - 1mm (20pc)
 €26.00 - 2mm (10pc)

Playsafe Triple Foils pack of 5 foils


€43.50 Single Colour
 €51.60 Two Colour
 €54.10 Three Colour
 €56.50 Four Colour

Playsafe Triple Foils/Pre Laminated Single Foils
 €14.75 Single Colour €17.20 Three Colour
 €16.40 Two Colour €17.98 Four Colour


Silensor-sl Kit €39.95



Silensor-sl Parts Card €15.00



Hands on course available



**One of the worlds leading snoring
appliances**

★ FDA Approved ★



Dentaclean denture cleaning agent



Dentaclean denture cleaning agent
1000 ml concentrate to obtain 11 liters of ready-to-use solution
REF 520 0099 2

Concentrate for easy removal of plaque, tartar and coatings on dentures. The liquid concentrate yields the 11fold quantity and therefore ensures high economic efficiency.

€58.00



Up to now it has always been very difficult to remove tartar from orthodontic screws and clasp edges.



Dentaclean denture cleaning agent allows for easy removal of tartar from orthodontic appliances without damaging filigree metal parts.

Before



After



Dentaclean impression and denture disinfectant



Disinfecting with Dentaclean impression and denture disinfectant avoids the transmission of viruses, bacteria and fungi from the patient to the laboratory. The concentrate yields 10 liters of ready-to-use solution which is highly effective and has a pleasant mild odour.

Dentaclean impression and denture disinfectant
1000 ml concentrate to obtain 10 liters ready-to-use solution including 25 shipping bags
REF 520 0100 6

€67.00



Pathogens can be transmitted to the laboratory with impressions.



After the use of Dentaclean impression disinfectant, active viruses, bacteria and fungi can no longer be detected.

Tested and approved by the Institute for clinical hygiene and infection control, Giessen.

Dentaclean mixing fluid for pumice powder



Protects against germs.

Dentaclean mixing fluid for pumice powder

- Remains moist and free of germs for two to three weeks without having to be remixed.
- Contains skin-care additives to protect employees' hands.
- Contains natural odours which still smell fresh after several weeks.
- Mixed polish adheres to the brush and restoration better so that the pumice splatters less. This saves time when polishing as the pumice slurry does not have to be applied repeatedly



In wet pumice powder disease microbes are present. Disinfection action occurs within one hour.



The nurturing effect for the skin is achieved by skin nurturing additive.

Dentaclean mixing fluid for pumice powder
5000 ml
REF 520 0099 8

Dentaclean mixing fluid for pumice powder
1000 ml
REF 520 0099 9

€24.50
1,000 ml

€98.00
5,000 ml

Application:

Simply mix the pumice powder with Dentaclean mixing fluid for pumice powder, do not add water. This way the pumice powder stays moist for two to three weeks.



MESTRA®

Automatic Alginate and Plaster Mixer

€595

With 6 FREE bags
of alginate



MESTRA®

Vibrating
Tables

Round, Rectangular
Square

R-080030

€225



R-080025

€245



Flexistone® Plus

Insulating and modelling
material, A-silicone based

€100.00 (5lt)

€25.00 (1lt)
**Separating
Fluid**



DETAX IDENTAL

€89.95 (2 x 160ml)



VIRTUAL TRAY ADHESIVE

10ml
Liquid

€9.95



**ivoclar
vivadent:**
passion vision innovation

Virtual

Vinyl Polysiloxane
impression material

- Excellent adaptation to moist oral tissues
- Outstanding physical properties
- Fresh mint scent
- Extra light body, light body monophase, heavy body
- Regular or Fast Set

**ivoclar
vivadent:**
passion vision innovation



CARTRIDGES

€31.95 (2 x 50ml)

TUBS

€56.95 (2 x 300ml)



SR Ivolen

Non-stick multi purpose
special tray material

€60.00 2 x 500g Powder

€26.00 500ml Liquid

KaVo. Dental Excellence.

Dental Laboratory Equipment FLEXspace

Customised laboratory furnishings

FLEXspace means individuality



The FLEXspace offers value

KaVo function and quality at attractive value for money prices. Economical for efficiently working in the long term. Modular design is a safe investment as when you move, the furnishings can be moved and set up in a new situation, and the table height is adjustable to suit new circumstances.

KaVo SENsit® Labor- und Bürostuhl / Laboratory- Office Chair

Stoffmuster / Fabric sample



Dust Extraction Systems

Single user, multi-user and central dust extraction systems for your dental laboratory and for your health.



KaVo MASTERspace

Perfect planning of the workstations is a pre-requisite of high efficiency. MASTERspace offers perfect quality, premium design, and optimal functionality.



**PLEASE CALL
FOR DETAILS**

Laboratory handpieces

Comfortable and powerful laboratory drives – including universal control unit.



ARCUSdigma supports therapy & diagnosis

ARCUSdigma, instrument-based functional diagnosis: a milestone towards the "digital office" and ideal support for your diagnosis and therapy. It was never that easy!



Dental Articulators KaVo PROTAR evo

The innovative evolution of dental articulators is called KaVo PROTARevo, which sets standards in the dental articulators market and dental diagnostics



**TASSO SAMBANIDIS 0872847795
ARDAGH DENTAL (01)8551600
info@ardaghdental.ie**



Classic Model Trimmer

R - 080082 RH 3000



€695

3000 rpm Silicon Carbide Disc Trimmer

Brushless Micro Motor



€650

Digital watch. Strong torque at low speed. Vertical / horizontal placement. RPM: MAX 50,000 rpm / Torque: 7.8 N. cm. Input: 100, 110, 210-230, 240 V (50/60 Hz) Output: 230 W. Right & Left turning ability. Hand / Foot speed control system. Self-diagnosis, Error display, Work memory storage function. Foot speed control



Water Bath

€295



High Capacity Wax Heater

€325



Curing Light Box



€295



€195

Induction Heater

TASSO 087-2847795

INFO@ARDAGHDENTAL.IE

ARDAGH DENTAL (01)8551600

ALL PRICES EXCLUDE VAT & OE

Zhermack
BEYOND INNOVATION

Gingifast
Laboratory A-Silicone for
gingival reproduction for
models in fixed
protheses



GINGIFAST ELASTIC OR RIGID
2 x 50ml cart, 12 mixing tips,
12 fine tips, separator

€41.00

Zhermack
BEYOND INNOVATION

Zetalabor

C-silicone laboratory putty
for duplication of denture
moulds



5kg - **€87.00**
10kg - **€162.00**



Indurent Gel
€12.99

Zhermack
BEYOND INNOVATION

Chromatic Alginate



€99.00

12 x 500gm

Zhermack
BEYOND INNOVATION



12 x 500gm

€75.75

Zhermack
BEYOND INNOVATION



12 x 500gm

€99.00

Elite LC Tray
Light-curing resin plates



€49.95



€64.95

Acrytemp

Bisacrylic self-curing composite
resin for temporary C&B
preparation