

**Ardagh**

Dental  
Supply  
Co. Ltd.




**BEST PRICES  
ON CAD BLOCKS**

# SOLUTIONS FOR DIGITAL DENTISTRY



**GUARANTEED  
DELIVERY**  
ON POPULAR SHADED  
e.max & Empress  
**CAD BLOCKS**



**SPECIAL OFFERS  
ON ALL CAD CAM  
PRODUCTS**

 **Roland**

 **bredent**

 **e.max**

 **3M Lava**  
Precision Solutions

 **Whip Mix**  
more reasons. one source

 **Komet**

 **Ivoclar  
vivadent**

**TASSO SAMBANIDIS 0872847795**  
[info@ardaghdental.ie](mailto:info@ardaghdental.ie)

**ARDAGH DENTAL (01)8551600**  
ALL PRICES EXCLUDE VAT & OE

## Restorative

# Make it IPS e.max®... because it works!


## IPS e.max® CAD

The leading CAD/CAM Lithium Disilicate Ceramic

- Veneers, inlays/onlays
- Crowns (monolithic or CAD-on technique, i.e. with ZrO<sub>2</sub> framework)
- 3-unit anterior and premolar bridges
- 3- to 4-unit bridges (CAD-on technique)
- Implant superstructures crowns (monolithic or CAD-on technique)
- Implant superstructures bridges (CAD-on technique)
- Hybrid abutments and hybrid abutment crowns



What would *dentists* put in *their* mouths?



# Make it IPS e.max®... because it matters!

IPS e.max CAD

### A14 – Innovative

Ideal for aesthetic hybrid abutments

The innovative characteristic of A14 blocks is based on the special interface, e.g. for Sirona TiBase. The individually fabricated hybrid abutments are permanently luted to titanium bases by means of Multilink® Hybrid Abutment.



IPS e.max CAD

### A16 – Efficient

Ideal for screwed-in hybrid abutment crowns

A16 blocks feature a special interface, e.g. for Sirona TiBase, and are particularly suitable for screwed-in hybrid abutment crowns. This renders the cumbersome step of removing excess luting material redundant. A16 blocks are available in the translucency level LT as well as in the nine well-established A–D shades and one Bleach shade.



IPS e.max CAD

### C16 – Large Classic

Ideal for large monolithic crowns

If C14 blocks are not sufficient, C16 blocks will do the trick. C16 blocks are used for large crowns (e.g. implant-supported crowns). They are available in the translucency level LT as well as in the nine well-established A–D shades and one Bleach shade.



IPS e.max CAD

### B32 – Long

Ideal for monolithic three-unit bridges

B32 blocks for three-unit bridges in the anterior and premolar region expand the treatment options within the digital processing chain. These blocks are available in the translucency level LT as well as in the nine well-established A–D shades and one Bleach shade.



## BLOCK PRICE LIST

**I12 & C14**

**€78.00 (PK / 5)**

**C16**

**€124.95 (PK / 5)**

**B32**

**€150.00 (PK / 3)**

**B40 & B40L**

**€169.00 (PK / 3)**

**A14**

**€150.00 (PK / 5)**

**A16**

**€185.00 (PK / 5)**



# Restorative



## IPS e.max® CAD

### The clinically proven solution for all your CAD/CAM needs

This unrivalled system combines highly aesthetic quality with exceptional user-friendliness.


- Comprehensive indication spectrum
- Scientifically and clinically proven
- Wide availability due to the CAD/CAM partner network (fabrication in-house and by authorised milling centres)
- Tried-and-tested cementation materials are coordinated with lithium disilicate (LS2)



#### NEW Blocks now available

IPS e.max CAD MT,  
IPS e.max LT A14 (S) & (L)

Buy any 2 packets of IPS e.max® CAD,  
receive an IPS e.max Mousepad FREE!



#### IPS e.max® CAD Crystall./Shades

Dentine stains in paste form that are optimally coordinated for use with IPS e.max CAD.



€39.95 (1x 3gm)

#### IPS e.max® CAD Crystall./Stains


Intensive stains in paste form for the reproduction of natural characteristics.



€18.95 (1x 1 gm)

#### IPS e.max® CAD Crystall./Glaze Paste

Special glaze paste which has been developed for use with fully anatomical IPS e.max CAD restorations.



€39.95 (1x 3 gm)

#### IPS e.max® CAD Crystall./Glaze Spray

Glaze spray to be used with fully anatomical IPS e.max CAD restorations.



€64.95 (1x 270ml)

## NEW Programat® CS2

### Your new glazing and crystallisation furnace

The second generation of the Programat CS series is the ideal ceramic and crystallisation furnace for dentists.

The Programat CS2 has been specially developed for the crystallisation of IPS e.max CAD restorations. It is easily operated by means of the new colour touch screen and the proven membrane-sealed keypad.

The Programat CS2 is equipped with 20 individual programmes as well as pre-defined Ivoclar Vivadent programmes, such as the Speed programme, which crystallises IPS e.max CAD restorations in less than 15 minutes.

- QTK2 muffle technology with SiC bottom reflector
- Easy operation; specially developed for use in the dental practice
- Automatic double-range temperature calibration (ATK2)
- Remote diagnostics via Internet, e-mail and data log files
- Software updates via USB flash drive
- OSD display
- Integrated tong holder
- Double valve vacuum technology



The furnace that opens up  
a whole new world.  
For dental practices.

**SPECIAL  
OFFER**

**CALL FOR DETAILS**



# Restorative



## IPS Empress® CAD

### Mill – Polish – Done!

- Twenty years of clinical experience with more than 40 million seated restorations
- Leucite glass-ceramic (160 MPa) for exceptional aesthetics
- Comprehensive shade and size range for every patient situation
- Choice of true-to-nature chameleon effect, lifelike brightness value or natural shade gradient
- Polishing with OptraFine – no glazing necessary

**CLINICALLY  
PROVEN**



### HT blocks

The IPS Empress CAD HT blocks with high translucency are mainly used for the milling of smaller restorations (e.g. inlays), due to their lifelike chameleon effect.

### LT blocks

Due to their chroma and brightness value, lower translucency (LT) blocks are ideally suitable for the fabrication of larger restorations (e.g. partial crowns, crowns).



### Multi blocks

The innovative polychromatic IPS Empress CAD Multi Blocks result in restorations which are distinguished for a true-to-nature shade and fluorescence transition from dentine to incisal, even without characterisation. They thus provide maximum aesthetics and true-to-nature appearance.

### IPS Empress® Universal Shades

Available in 15 Chromascop and 9 A-D shades.

**€23.00** (1 x 3gm)



### IPS Empress® Universal Stains

Available in 14 different shades and used to imitate natural characteristics.

**€11.00** (1 x 1gm)



### IPS Empress® Universal Glaze

Glaze paste to glaze all IPS Empress CAD restorations.

**Glaze & Stain Liquid**

**€15.00** (1 x 15ml)



**Glaze Paste  
€38.00**  
(1 x 3gm)

## Optrafine

**KIT - €99.00**

The OptraFine polishing instruments and high-gloss polishing paste feature a higher than average concentration of diamond particles. This ensures exceptional smoothness and unsurpassed aesthetic lustre of the ceramic surfaces. You can therefore achieve perfect polishing results in no time and without compromise.

- Unlike other polishing systems, the high-quality OptraFine polishers can be re-used up to 10 times.



## BLOCK PRICE LIST

**HT / LT I10 / I12 / C14 - €49.00 (PK 5)**

**MULTI C14 - €65.00 (PK 5)**

**MULTI I12 - €65.00 (PK 5)**

Blocks include, IPS e.max, Telio CAD and IPS Empress CAD. Cheapest block packs will be FREE. Offer only available when all placed as one order.

### PRODUCT HIGHLIGHT

A14 S & L CAD Blocks now available in LT: A1, A2, A3, A3.5, B1, B2, C1, C2 & D2



**TASSO SAMBANIDIS 0872847795**

[info@ardaghdental.ie](mailto:info@ardaghdental.ie)

**ARDAGH DENTAL (01)8551600**

ALL PRICES EXCLUDE VAT & OE



## Cementation Solutions

### Multilink Automix



**Multilink Hybrid Abutment Cement**  
**Refill €129.95**



### Variolink Esthetic



**DC System Kit - €320.00**  
**DC Promo Kit - €240.00**



**System Kit**

**€184.00**

**Refills**

**€129.95**

**Bottles A & B**

**€58.00**



### Telio ®CS Link



**€72.00**

Dual-curing, translucent, temporary composite cement for the esthetic and provisional, cementation of temporary restorations

### Telio CS Desensitizer

Telio CS Desensitizer is a clear solution for reducing and preventing dentinal hypersensitivity and postoperative sensitivity.



**€48.95** (5gm)

### IPS Ceramic Etching Gel

**€25.00** (5gms)




**Monobond Etch & Prime**  
Self-etching glass-ceramic primer  
**World's first self-etching glass-ceramic primer shortens the pretreatment of dental restoration**

**€89.95**

### SpeedCEM Plus

SpeedCEM is a self-adhesive, self-curing resin cement with light-curing option.



**€99.95** (9gms)

### Ivoclean

Proper cleaning after try-in

**€20.00** (5gms)




# ivoclar vivadent®

passion vision innovation

## Accessories

### CAD CAM SOLUTIONS



**OptraGate Refill**  
**€84.95**



**IPS e.max CAD Crystallization Tray**  
**€99.00**

**IPS Contrast Spray  
Chairside (50ml)**  
**€42.95**



**IPS Object Fix Flow & Putty**  
**€16.95** (10gm)



**OptraStick** **€18.95** (50 Sticks)



**"ATK 2"**  
**Calibration Sticks**  
**€180.00**



**Virtual CADbite registration**  
**€38.95** (2 x 50ml)



**Ceramic Repair Kit** **€195.00**



**Natural Die Material Kit** **€450**  
**Refill** **€39.95**



**IPS Contrast  
Spray Labside**  
**€89.95** (270ml)



**Red Marten Modelling Brushes Kit**  
**€89.95**

**Refill Sizes**  
00 - **€8.95**  
0 - **€8.95**  
1 - **€8.95**  
2 - **€8.95**  
4 - **€8.95**  
6 - **€17.50**  
10 - **€45.80**  
G2 - **€9.50**  
G4 - **€11.65**





# CAD CAM SOLUTIONS



**€89.95** (Pk 5) I12 & C14 HT/LT  
SHADES A1,A2,A3,A3.5,B1,C2,D2,BLEACH



**€495.00** Cojet Sand Blaster



**€185.00** RelyX Unicem  
Aplicap (50 Caps)



**ASK US!**



**RelyX Ultimate Cement**

**Trial Kit €149.00**

**Refill €135.00**

**Scotchbond Universal €64.95**

**Sinfony Shades**

A micro-hybrid, indirect laboratory composite



**Astringent Retraction Paste**

**€55.00** (25 Caps)

**€189.00** (100 Caps)



**Sinfony Magic Refill**

**€62.95** Per Refill




**TASSO SAMBANIDIS 0872847795**  
[info@ardagh dental.ie](mailto:info@ardagh dental.ie)

**ARDAGH DENTAL (01)8551600**


ALL PRICES EXCLUDE VAT & OE

**bredent**  
More than 40 years of dental innovations



**BioHPP**

BioHPP® is a PEEK-based, ceramic-reinforced high performance polymer.



## breCAM BioHPP



### High Performance Peek Material

Available in Ø 12 mm -16 mm -20 mm -24 mm

breCAM BioHPP €275.00 (12mm)

breCAM Bio HPP €320.00 (16mm)

breCAM Bio HPP €340.00 (20mm)

breCAM Bio HPP €355.00 (24mm)

## breCAM Resin



### High Performance PMMA

Available in Ø 16 mm -20 mm -24 mm

breCAM.resin A Shade €56.50 (16mm)

breCAM.resin B Shade €64.95 (20mm)

breCAM.resin C Shade €69.95 (24mm)

breCAM.resin transparent €59.95 (20mm)


## breCAM Wax



### breCAM. wax

Ø 98,4 x 20 mm €44.95 (2 piece)

## breCAM multiCOM



### High Performance PMMA with ceramic fillers (wet/dry mill)

Shades A1 / A2 / A3 / A3.5 / B2

16mm €79.95 / 20mm €89.95




### Worlds first individual Physiological Hybridabutment

3 Pieces €230.00

6 Pieces €395.00

9 Pieces €585.00

12 Pieces €730.00



**bredent**

**Diagen-Turbo-Grinder**

**Assortment €109.00**

**Refills Fine 2-pcs/pkg**

**€32.95**

**Refills Coarse 2-pcs/pkg**

**€34.95**



### breCAM HIPC "High Impact Polymer Composite"

Shades Bleach A1, A2, A3, A3.5, A4, B2, B3, C2, C3, D2, D3, trans.



**breCAM.HIPC**

**TASSO SAMBANIDIS 0872847795**


**ARDAGH DENTAL (01)8551600**

**info@ardaghdental.ie**

**ALL PRICES EXCLUDE VAT & OE**



4637.000




€45.95

★ Set for minor corrections and polishing of ZrO<sub>2</sub>

● ○	ZR379L.315.014	1	
● ○	ZR8379L.315.014	1	
● ○	ZR850.314.016	1	
● ○	ZR8850.314.016	1	
	94020C.204.040	1	
	94020F.204.040	1	
	94022C.204.060	1	
	94022F.204.060	1	

## Ceramic Veneers | Kit 4388



€39.50

★ Ceramic Veneers - a kit of coordinated diamond instruments for the conservative, controlled preparation of ceramic veneers

Depth reduction	Shaping
868B.314.018	868.314.012
868B.314.020	868.314.016
	379.314.023
Fine correction	Finishing
852EF.314.014	8868.314.012
379EF.314.023	8868.314.016
	8379.314.023

## 4562 | 4562ST



€48.95


★ Expert Set for ceramic inlays and partial crowns  
Set 4562ST  
In a bur block suitable for sterilisation

## Content of Set 4562/4562ST

6847KRD.314.016	8847KR.314.016
959KRD.314.018	8959KR.314.018
845KRD.314.025	8845KR.314.025
8862.314.012	8862.314.016
801.314.023	8379.314.023

€37.95

## Guide Pin | Kit 4410



★ for all-ceramic preparations with guide pin instruments

8845KR.314.018	8372P.314.023
379.314.023	8372PL.314.023
H50A.314.010	8847KR.314.023
9653.204.060	8845KR.313.025
6830LGK.204.014	8846KR.314.018

## 4573 | 4573S



€47.95

★ Expert set for ceramic crown preparation

## Contents of sets 4573S/4573

6837KR.314.012	6836KR.314.014
6856.314.021	8856.314.018
6856.314.018	8856.314.021
6856.314.012	
6379.314.023	8379.314.023

## Trim Gingival Tissue | CeraTip™



€95

Set 4561

## Crown preparation | Set 4384A



€79.95

★ Diamond instruments with guide pin for controlled tissue-preserving crown preparation

878KP.314.018 ○ <sub>rev</sub> 40.000 rpm	878KP.314.021 ○ <sub>rev</sub> 40.000 rpm	856P.314.018 ○ <sub>rev</sub> 40.000 rpm	856P.314.021 ○ <sub>rev</sub> 40.000 rpm
8878KP.314.018 ○ <sub>rev</sub> 20.000 rpm	8878KP.314.021 ○ <sub>rev</sub> 20.000 rpm	8856P.314.018 ○ <sub>rev</sub> 20.000 rpm	8856P.314.021 ○ <sub>rev</sub> 20.000 rpm
68878K.314.012	68856.314.018	6379.314.023	
68878K.314.016	8856.314.018	661.314.420	



The CeraTip is suitable for use in various sectors of mucosa surgery, such as:

- Exposing intra-osseous implant sites ( g. 1)
- Dilatation of the sulcus following a crown core preparation ( g. 2)
- Exposing deep cavities in the neck of the tooth
- Exposing impacted teeth
- Papillectomy / Removal of hyperplastic gingiva



New **DWX-51D** Dental Mill

Two New Dental Mills  
Designed to Meet Your Needs




New **DWX-4W** Wet Dental Mill

## **NEW** The DWX-4W Wet Dental Mill

Introducing the DWX-4W Wet Dental Mill, the latest in Roland DG technology that goes beyond dry milling to quickly and precisely cut glass ceramics and composite resins. The DWX-4W produces crowns, bridges, inlays and onlays, and veneers at high-speed with its Jaeger DentaDrive spindle and simultaneous 4-axis milling and offers unattended milling convenience and increased efficiency with the option of loading up to three standard pin-type blocks at a time. Featuring Roland DG's famous reliability and product support, it's the ultimate in affordable and dependable milling technology.

**Call for Details**



The Latest in DWX Innovation  
**DWX-4W**

**EASY  
SHAPE**  
Roland  
DENTAL SOLUTION

## **NEW** The DWX-51D 5-Axis Dental Milling Machine

The DWX-51D 5-Axis Dental Milling Machine is designed for effortless production of dental prosthetics. Built for precision, error-free milling, you don't need any machine expertise to operate it. With touch-button controls and automated milling features, the production of high-quality bridges, crowns, inlays, onlays, abutments and other restorations out of 98 mm discs and blocks has been made virtually foolproof.



**TASSO SAMBANIDIS 0872847795**

**ARDAGH DENTAL (01)8551600**

[info@ardaghdental.ie](mailto:info@ardaghdental.ie)

**ALL PRICES EXCLUDE VAT & OE**





**Ardagh** | Dental  
Supply

*more reasons. one source*

**Whip Mix Veribuild 100**  
Ideal for printing patterns, temporaries  
and flexible gingiva masks

**3D**  
**PRINTERS**

**€5,995**



**Whip Mix Veribuild 200**

- Your choice for Restorative applications
- Very versatile
- Ideal for printing patterns, drill guides, models, implant trays, temporaries and flexible gingiva masks

**€23,950**



**CALL FOR DETAILS**

**Whip Mix Veribuild 300**

- Your choice for Ortho and all crown and bridge applications
- Largest build plate
- Ideal for printing patterns, drill guides, models, implant trays temporaries and flexible gingiva masks

**€29,995**



**TASSO SAMBANIDIS 0872847795**  
**ARDAGH DENTAL (01)8551600**

# Chairside Scan, Design, Mill

A truly open system delivered by truly open partners



## SCAN

With 3Shape TRIOS3

## DESIGN

With 3Shape Practice Lab

## MILL

With Roland DWX-4W

### Chairside Mill with TRIOS3 CART

3Shape TRIOS3 Cart  
HP Elite 800G1 PC + Monitor  
3Shape Practice Lab Software  
Roland DWX-4W wet Mill  
Sum3D millbox CAM software  
IPS Programat CS2 Furnace (w/vacuum pump)  
Starter Kit (burs, blocks, crystallisation kit)  
Delivery, Installation & Training  
1yr warranty

RRP: €70,663

Introductory Offer: €62,050

### Chairside Mill with TRIOS3 POD

3Shape TRIOS3 POD  
HP Omen 15 Laptop  
HP Elite 800G1 PC + Monitor  
3Shape Practice Lab Software  
Roland DWX-4W wet Mill  
Sum3D millbox CAM software  
IPS Programat CS2 Furnace (w/vacuum pump)  
Starter Kit (burs, blocks, crystallisation kit)  
Delivery, Installation & Training  
1yr warranty

RRP: €66,863

Introductory Offer: €59,850

### Chairside Mill with TRIOS Mono POD

3Shape TRIOS Mono POD  
HP Omen 15 Laptop  
HP Elite 800G1 PC + Monitor  
3Shape Practice Lab Software  
Roland DWX-4W wet Mill  
Sum3D millbox CAM software  
IPS Programat CS2 Furnace (w/vacuum pump)  
Starter Kit (burs, blocks, crystallisation kit)  
Delivery, Installation & Training  
1yr warranty

RRP: €55,863

Introductory Offer: €50,150

Contact ESM Digital Solutions at (01) 808 4446 or Ardagh Dental Supplies at (01) 855 1600 for further information and to arrange a presentation at your practice.

- 1) the customer is responsible for ensuring that compressed air is provided according to specification (continuous supply of dry, oil free air, min 6bar @114 l/min)
- 2) Prices do not include VAT or any other taxes and duties.
- 3) Add-ons priced separately.
- 4) Training is provided at the customers site and is included in the above prices.
- 5) An annual subscription fee is payable and ensures on-going software licences, upgrades and remote support and training.
- 6) Standard warranty terms 1yr. Extended warranty options are available.
- 7) Introductory prices available for a limited time only.

**Ardagh** Dental Supply

**Roland**

**ESM**  
DIGITAL SOLUTIONS

**3shape**

**EUROCAST**